

## The US Pro 2000™ 2nd Edition

The **US Pro 2000™** 2nd Edition is a clinical grade therapeutic ultrasound device that generates deep heat through ultrasound energy. The pulsed sound waves travel deep into the tissue and induce vasodilation, which helps increase blood flow to the treated area and is found to help relieve pain and reduce muscle spasms. Equipped with a head warming feature on an ergonomically designed sound head for increased patient comfort.

With three output intensities and three treatment time selections, the US Pro 2000 2nd Edition can provide a wide range of therapy treatments. Sound head is watertight for use in underwater therapy (for 30 mins at depth of 1 meter). Complete with hard carrying case, ultrasound gel, AC/DC adapter, instruction manual and a quick start guide.

Part# 710552-001

### FEATURES

- Operates in continuous mode (100% Duty Cycle).
- Head warming feature pre-heats up to 95°F the sound head applicator for increased patient comfort.
- Sound head only (aluminum portion – IPX7) can be used for sub-aqua therapy in water no more than 1 meter deep.
- Three output intensities selections.
- Three treatment times (5 min, 10 min or 15 min).
- Small, ergonomic design that is easy to hold and portable.
- Comes complete with device, ultrasound gel, AC/DC adapter, hard carrying case, quick start guide and instruction manual.
- UPC Code: 0 92237611555
- One year warranty

Distributed by:  
**TruMedical Solutions, LLC.**  
Phone: 877.882.7844  
www.tru-medical.com

© 2016 TruMedical Solutions, LLC.

