



## TREADMILL T635M

SPORTSART REHAB TREADMILLS ARE DESIGNED SPECIFICALLY FOR CARDIAC AND PHYSICAL THERAPY MARKETS AND COME LOADED WITH UNIQUE FEATURES, PROGRAMS AND USER AMENITIES TO HELP TREAT MORE PATIENTS MORE EFFECTIVELY.

- EQUIPPED WITH MEDICAL/REHAB SPECIFIC FEATURES SUCH AS REVERSE SPEED UP TO 3 MPH AND EXTENDED HANDRAILS.
- MYFLEX PLUS™ DECK CUSHIONING PROVIDES 30% MORE SHOCK ABSORPTION AND CARRIES A LIFETIME WARRANTY.
- MEETS UL REQUIREMENTS FOR CURRENT LEAKAGE IN MEDICAL SETTINGS.
- LARGE 61 X 22 IN BELT PROVIDES AMPLE RUNNING SPACE AND IS CONSTRUCTED FROM CARBON-WEAVE, LOW FRICTION MATERIAL.

OVERALL WEIGHT	DIMENSIONS	STEP UP HEIGHT	RUNNING SURFACE	INCLINE RANGE	SPEED RANGE	MAXIMUM USER WEIGHT	DRIVE MOTOR
353 LBS	84 x 39 x 57IN	7.5"	61 x 22IN	0 TO 15%	-3 TO 10 MPH	400 LBS	3.2 HP DC

### FEATURES AND SPECIFICATIONS

#### FEEDBACK.

CALORIES , SPEED , TIME , DISTANCE , CAL/HR , METs , PACE , INCLINE , HEART RATE , WEIGHT LOSS ZONE , CARDIO ZONE , ACTIVE ZONE

#### WORKOUT PROGRAM

QUICK START , MANUAL , RANDOM (X10,000) , INTERVAL (X3) , GLUTE (X2) , HILL (X3) , 5K/10K , FAT BURN , HRC (X3) , FITNESS TEST (BRUCE, GERKIN)

#### POWER REQUIREMENT

15A/120V/60Hz - DEDICATED CIRCUIT

#### ACCESSORIES

- CSAFE PORT
- MEDICAL HANDRAILS
- OPTIONAL WIRED CONTROLLER
- FORWARD AND REVERSE SPEED (03 TO 10MPH)
- O.I SAFESTART™ SPEED (FORWARD & REVERSE)
- CONTACT AND TELEMETRY HR
- MULTI-PLY 61 X 22IN BELT
- USB PORT (X2) FOR DEVICE CHARGING
- OPTIONAL 3-SPEED FAN
- ECO-GLIDE™ AUTO LUBRICATION SYSTEM
- CENTER HANDLE BAR
- OPTIONAL SA WELL+™
- DUAL HEADPHONE JACK WITH BUILT IN VOLUME CONTROL

#### WARRANTY

##### LIGHT COMMERCIAL

- FRAME: LIFETIME
- DRIVE MOTOR: LIFETIME
- PARTS: 5 YEARS
- LABOR: 3 YEARS