



2020 TRU™ BRAND PRODUCT CATALOG



Partnering With The Provider
To Meet The Needs Of The Patient

ORDERING INFORMATION

Order Options:

Order by phone, fax, mail
or web.

Phone Toll Free:

1-877-882-7844

Hour of Operation:

Monday thru Friday
8:30 am to 6:00 pm EST.

Fax Toll Free:

1-877-882-7845


24 hours a day, 7 days a week.

Mail To:

TruMedical Solutions, LLC.
PO Box 1869
Collegedale, TN 37315

Order By Web:

www.tru-medical.com



We Accept All Major Credit Cards:



TABLE OF CONTENTS

Hot & Cold.....	2-3	Training Stairs.....	11-12
Electrodes, Lotion & Gel.....	4	Tables, Parallel Bars, Neck Pro	13-16
Band, Pulleys & Tubing.....	5-6	Taping, Work & Mat Table.....	16-19
Weight Racks, Stools, Mirror.....	7-9	Bolsters & Wedges	19-20
Therapy Putty	9	Pillows, Fasteners, TruBIA.....	20-21
Pillows, Rolls, Cushions	10-11, 25	Stockinette, Tubular Bandage	22
Kiosk, Hi-Low Stools.....	11	Cervical Collars, Belts, Dynamometer.....	23-25

Helpful Information At A Glance

Look for these icons throughout the catalog to easily identify special product features...



Contains No Latex



Ships By Freight Carrier



Requires Some
Assembly

Signature
Series

Signature Series
Silkscreen



ORDER ONLINE AT WWW.TRU-MEDICAL.COM



WEBSITE FEATURES

- Always open to see your pricing
- Place and track orders
- See your quotes
- Print a statement
- View all your transactions, and much more...
- We store all your information to make it simple to shop at
- We offer 9 different categories of product to meet all your clinic needs.

VISIT US TODAY AT WWW.TRU-MEDICAL.COM



TruBIA™ Body Composition Analyzer

A great new product that measures impedance for body composition using surface electrodes. Offers accurate and repeatable results that will articulate total water analyzes, hydration and wasting analysis. Competitively priced.

770177-001	TruBIA™ Body Composition Analyzer.....	\$2,195.00
770177-002	Electrodes, 100/Package.....	19.95
770177-003	Lead-Wire, Pair	70.00

More Information on TruBIA™ Pg. 21

Hot Packs



710103-001	Standard 10" x 12"	\$12.25
710103-002	Cervical 24" Long	13.95
710103-003	Oversize 15" x 24"	25.50
710103-004	Knee & Shoulder 10" x 20"	24.50
710103-005	10" x 18"	20.50
710103-006	10" x 24"	23.50
710103-008	Half 5" x 12"	9.25

Tru™ Hot Packs

- Up to 30 minutes of deep moist heat
- Blue tabs on Hot Packs for easy rotation in busy clinics
- Each pack provides hundreds of treatments
- Can be used in the clinic and at home
- Great for treatment of arthritis, strained ligaments & muscle spasms
- Variety of sizes accommodate different areas of the body
- 12 month warranty
- Private Label Terry Covers Have Minimum Order 36 Mix & Match!

710103-011	"Private Label" Hot Pack, Standard 10"x 12"	\$12.25
710103-022	"Private Label" Hot Pack, Cervical 24" Long	13.95
710103-033	"Private Label" Hot Pack, Oversize 15"x24"	25.50
710103-044	"Private Label" Hot Pack, Knee & Shoulder 10" x 20"	24.50
710103-055	"Private Label" Hot Pack, 10"x18"	20.50
710103-066	"Private Label" Hot Pack, 10"x24"	23.50
710103-088	"Private Label" Hot Pack, Half 5" x 12"	9.25
710103-100	"Private Label" Hot Pack, Setup Charge.....	80.00

Tru™ Terry Covers

- Designed to cover hot packs to insulate and protect the patient from excessive heat exposure.
- Available in durable, long-lasting all-terry and quick-drying foam-filled styles
- For high institutional use
- Hook and loop closures prevent unwrapping and slipping
- Wash and drip dry for best results
- 12 month warranty
- Private Label Terry Covers Have Minimum Order 36 Mix & Match!

All Terry

710101-001	Standard.....	\$24.50
710101-002	Cervical	19.95
710101-003	Oversize.....	33.25
710101-011	"Private Label" Standard All Terry Cover.....	24.50
710101-022	"Private Label" Cervical All Terry Cover	19.95
710101-033	"Private Label" Oversize All Terry Cover.....	33.25



ⓐ All Terry



ⓑ Foam Filled

Foam Filled

710102-001	Standard.....	\$18.95
710102-002	Cervical	18.50
710102-003	Oversize.....	26.50
710102-011	"Private Label" Standard Foam Filled Terry Cover	18.95
710102-022	"Private Label" Cervical Foam Filled Terry Cover.....	18.50
710102-033	"Private Label" Oversize Foam Filled Terry Cover.....	26.50

Cold Packs

TruTM *Signature Series*

Customize TruTM Brand Products With Your Clinic Name And Logo


TruTM Deluxe Cold Packs

- Provides soothing therapy for acute pain due to injury or inflammation
- Reusable design offers up to 30 minutes of cold therapy
- Conforms to any body contour
- Made of durable polyurethane (rubber) material for longer life
- Flexible filler will not freeze solid
- Ask customer service about our private label program
- 12 month warranty
- Made in the USA

Signature Series

A	710126-001	Standard, 10" x 13½"	\$18.95
B	710126-002	Cervical, 21" Long	18.95
C	710126-003	Oversize, 12½" x 18½"	30.50
D	710126-004	Half Size, 6½" x 11"	11.25
	710126-005	Quarter Size, 5½" x 6½" (Picture Not Shown)	9.15

A minimum of 36 Tru "Signature Series" Cold Packs are required. You may mix and match sizes to meet the minimum quantity. There is a one-time Die Charge of \$80.00 (710126-100) for the customized name and logo. Logo dimension size is 2" x 4".



710126-011	Private Label Standard, 10" x 13½"	\$18.95
710126-022	Private Label Cervical, 21" Long	18.95
710126-033	Private Label Oversize, 12½" x 18½"	30.50
710126-044	Private Label Half Size, 6½" x 11"	11.25
710126-055	Private Label Quarter Size, 5½" x 6½"	9.15
710126-100	Private Label Cold Pack, Die Charge	80.00



TruTM Clinical Cold Compression Kits

- Provides excellent cold plus compression therapy for up to 30 minutes
- Kit offered with size appropriate cold pack, compression wrap and straps
- Designed to provide evenly compressed adjustable cold therapy
- Clinical wrap is made of durable and cleanable coated material for clinic use

- Standard Kit treats shoulder, back, hip & knee
- Cervical Kit treats the neck
- Half-Size Kit treats elbow, wrist & ankle
- If you already own the cold packs, we offer the wraps & straps only to make your own kit
- Made in the USA

Cold Pack & Compression Wrap Kit

710471-001	Standard	\$47.90
710471-002	Cervical	47.90
710471-004	Half Size	34.50
710471-005	(2) Standard	71.40
710471-011	Standard w/X-Long Straps	52.50
710471-055	(2) Standard w/X-Long Straps	76.00

Compression Wrap Only

710472-001	Standard Cold Pack	\$28.95
710472-002	Cervical Cold Pack	28.95
710472-004	Half Size Cold Pack	23.25
710472-011	w/X-Long Straps for Standard Cold Pack	33.50



TruTM

TruTM Cold Pack Freezer Chilling Unit

- Includes 6 cold packs (oversize cold pack = (2) standard packs)
- Offered with 2 shelves to ensure even cooling of cold packs
- Dimensions: 22"D x 20.5"W x 32.5"H
- Adjustable temperature control
- Manual defrost function
- Does not include ice massage cups

710624-001	(3) Standard & (3) Cervical Cold Packs	\$550.00
710624-002	(2) Standard, (2) Cervical & (1) Oversize Cold Packs	550.00

710624-003	(6) Standard Cold Packs	\$550.00
710624-004	(2) Oversize & (2) Cervical Cold Packs	550.00

Electrodes




Self-Adhesive Reusable Electrodes

- Dual Layer Gel: increases number of applications
- High Quality Carbon: for uniform current distribution
- For All Electrotherapy Applications: TENS, NMES, FES, IF
- Made in the USA: ensures superior quality



710280-002
2" x 2" Square Cloth (4/pkg)

PO Box 1869 • Collegedale, TN 37315
1-877-882-7844 • www.tru-medical.com



Tru™ Self-Adhesive Reusable Electrodes

- Dual layer gel, increases number of applications
- Oversized border to facilitate proper and easy removal.
- High quality carbon film, for uniform current distribution
- For all electrotherapy applications: TENS, EMS, FES, IF
- Made the USA, ensures superior quality
- 4 electrodes per package

Cloth-Backed

710280-001	2" Round 4/Pkg		\$6.25
710280-002	2" x 2" Square 4/Pkg.		6.25
710280-003	3" Round 4/Pkg.		9.50
710280-004	2" x 3 1/2" Rect. 4/Pkg.		9.50


Foam-Backed

710281-001	2" Round 4/Pkg.		\$6.25
710281-002	2" x 2" Square 4/Pkg.		6.25
710281-003	3" Round 4/Pkg.		9.50
710281-004	2" x 3 1/2" Rect. 4/Pkg.		9.50




Self-Adhesive Reusable Cloth Electrodes

- Dual Layer Gel, increases number of applications
- High Quality Carbon, for uniform current distribution
- For All Electrotherapy Applications: TENS, EMS, FES, IF
- Zip lock bag for increased electrode life & storage



710582-004
2" x 3.5" Rectangle
4 PER PACKAGE

PO Box 1869 • Collegedale, TN 37315
1-877-882-7844 • www.tru-medical.com



MediFlex™ Reusable Electrodes

- Dual layer gel, increases number of applications
- High quality carbon film, for uniform current distribution
- For all electrotherapy applications: TENS, EMS, FES, IF
- Superior quality
- 4 electrodes per package

Cloth-Backed

710582-001	2" Round 4/Pkg.		\$5.95
710582-002	2" x 2" Square 4/Pkg.		5.95
710582-003	2.75" Round 4/Pkg.		8.95
710582-004	2" x 3 1/2" Rectangle 4/Pkg.		8.95



Tru™ Ultrasound Gel

- Excellent conductor for medical ultrasound treatments
- Will not irritate skin, stain clothes or damage transducers
- Water soluble, salt and alcohol free
- Offered in 5 liter container with one 8 oz. dispenser bottle
- Unscented

710065-001	5 Liter.....	\$16.50
710072-001	8 oz. Dispenser Bottle (empty).....	1.50



Tru™ Ultrasound Lotion

- A great multi-purpose ultrasound lotion
- Can also be used for massage therapy
- Water soluble, salt and alcohol free
- Offered in 5 liter container with one 8 oz. dispenser bottle
- Unscented

710066-001	5 Liter Lotion	\$17.50
710072-001	8 oz. Dispenser Bottle (empty).....	1.50

Exercise Band



Great Quality, Contains Latex



Tru™ Exercise Band

- Excellent progressive exercise system for improving strength, mobility and flexibility
- Fun colors for easy identification of five resistance levels
- Same great quality as higher priced name brands, yet easy on your budget
- Extremely versatile for unlimited exercise techniques and uses
- Offered in 50 yard dispenser box, 5ft. Length 10-Pack and Light/Heavy Packs
- Contains latex

50 yd. Resistive Exercise Band

720005-046	Level 1 Peach	\$46.95
720005-047	Level 2 Orange.....	56.95
720005-048	Level 3 Lime Green.....	66.95
720005-049	Level 4 Blueberry	73.95
720005-051	Level 5 Plum	77.95

5 ft. Exercise Band Dispenser 10-Pack

720485-046	Level 1 Peach	\$24.50
720485-047	Level 2 Orange.....	26.50
720485-048	Level 3 Lime Green.....	28.95
720485-049	Level 4 Blueberry	31.95
720485-051	Level 5 Plum	36.95

Exercise Band Pack

720486-001	(Includes One Each Level 1, 2 & 3)	\$8.25
720486-002	(Includes One Each Level 3, 4 & 5)	9.95



Versatile For Unlimited Exercise

Tru™ Latex-Free Exercise Band

- Excellent progressive exercise system for improving strength, mobility and flexibility
- Fun colors for easy identification of five resistance levels
- Same great quality as higher priced name brands, yet easy on your budget
- Extremely versatile for unlimited exercise techniques and uses
- Offered in 50 yard dispenser box
- Latex-Free

720417-046	Level 1 Peach	\$52.95
720417-047	Level 2 Orange.....	57.95
720417-048	Level 3 Lime Green.....	66.95
720417-049	Level 4 Blueberry	71.95
720417-051	Level 5 Plum	75.95



Tru™ Vertical Band Rack

- The Tru™ Vertical Band Rack dispenses (5) bulk rolls of exercise band
- Hand-crafted with high quality oak plywood
- All surfaces are finished with durable clear polyurethane
- Each roll easily mounts on a dowel rod to aid in dispensing band
- A storage shelf is offered for other supplies
- Dimensions: 43"H x 8"W x 8"D
- Made in the USA

720795-001 Band Rack w/Storage (holds 5 rolls) \$245.00



Great Quality Products
For Every Need.

Tru™ Latex-Free Exercise Tubing

- Excellent progressive exercise system for improving strength, mobility and flexibility
- Characteristics and performance of latex tubing, with the safety of non-latex allergy issues
- No latex, no powder, no odor
- Fun colors for easy identification of five resistance levels
- Extremely versatile for unlimited exercise techniques and uses
- Offered in 100 ft. and 25 ft. lengths
- Made in the USA



**LATEX
FREE**

720711-046	100 ft. Level 1 Peach.....	\$25.75
720711-047	100 ft. Level 2 Orange.....	42.25
720711-048	100 ft. Level 3 Lime Green.....	56.25
720711-049	100 ft. Level 4 Blueberry.....	72.25
720711-051	100 ft. Level 5 Plum.....	79.50
720719-046	25 ft. Level 1 Peach.....	9.25
720719-047	25 ft. Level 2 Orange.....	14.50
720719-048	25 ft. Level 3 Lime Green.....	18.50
720719-049	25 ft. Level 4 Blueberry.....	22.75
720719-051	25 ft. Level 5 Plum.....	26.50

Stretch Strap



Tru™



Tru™ Stretch Strap

- Stretching Aid for Triceps, Arms, Shoulders, Outer Thigh and Hips, Hamstrings, Quad, Calf, Low Back.
- Firm non-stretch webbing on one side and elastic stretchable loops on the other.
- Six feet of durable nylon webbing is complemented with tough elastic loops stitched in that provide anchoring position for feet and/or hands.
- Made in the USA

720640-001	Tru Stretch Strap.....	\$17.50
720640-002	Tru XL Stretch Strap.....	18.50

Pulleys



Tru™

Hand-crafted
with high quality
solid oak

Tru™ Deluxe Wall-Mount Pulley

- The Tru™ Deluxe Wall Pulley offers range of motion exercises for the shoulder in the clinic
- Hand-crafted with high quality solid oak
- A great transition to home shoulder pulley systems
- Made from 3/4" oak plywood, rope & pulley
- Dimensions: 18"H x 5"W x 20"D
- Installation hardware not included
- Made in the USA

720806-001	Tru™ Deluxe Wall Mount Pulley.....	\$110.00
------------	------------------------------------	----------



Deluxe
Pulley



Standard
Pulley



Buckle
Anchor



Non-scratch
Anchor

Tru™ Pulley Systems

- Cost effective devices designed to assist in shoulder rehabilitation
- Nylon webbing with buckle door anchor easily attaches to most conventional doors
- Non-scratch door model features a metal door anchor with foam pads
- The Tru™ Deluxe Pulley features a nylon pulley which glides smoothly, providing a more comfortable shoulder exercise
- Made in the USA

Standard Tru Pulley

720050-001	with Buckle Anchor.....	\$8.50
720050-002	with Non-Scratch Anchor.....	14.25

Deluxe Tru Pulley

720051-001	with Buckle Anchor.....	\$11.95
720051-002	with Non-Scratch Anchor.....	17.25

Weight Racks



Tru™ Multi-Purpose Weight Rack w/Mirror

- Beautiful craftsmanship on this multi-purpose weight rack
- White coated pegboard on front of rack for (9) therapy weights, includes (11) Pegboard safety hooks for therapy weights
- Includes (10) dumbbell hooks on left side and (10) hooks on the right side to store dumbbells
- Offers full view distortion-free glass mirror
- Two large shelves available for band/supplies storage
- Made in durable black or natural oak finish
- Offered with (4) 3" diameter easy-rolling casters
- Weights and bands sold separately
- Made in the USA

720587-001	Natural Oak.....	\$810.00
720589-001	Black	810.00
720750-001	Hooks for Weight Carts, 12 Pack	11.95



Tru™ Combination Weight Rack

- Beautiful craftsmanship on this combination weight rack
- White coated pegboard on front of rack for (9) therapy weights, includes (11) Pegboard safety hooks for therapy weights
- Includes (20) dumbbell hooks on back side to store dumbbells
- Two large shelves available for band/supplies storage
- Made in durable black or natural oak finish
- Offered with (4) 3" diameter easy-rolling casters
- Weights and bands sold separately
- Made in the USA

720586-001	Natural Oak.....	\$685.00
720588-001	Black	685.00
720750-001	Hooks for Weight Carts, 12 Pack	11.95



Tru™ Combo "Mini" Storage/Weight Rack

- Beautiful craftsmanship on this Combo "Mini" storage/weight rack
- White coated pegboard on front of rack for (9) therapy weights, includes (11) Pegboard safety hooks for therapy weights
- Includes (5) dumbbell hooks on left side and (5) hooks on the right side to store dumbbells, can hold (1) 1 to 10 lb. or (2) 1 to 5 lb. dumbbell sets
- The back of the rack offers four large storage shelves enclosed by beautiful oak doors to store all your clinic supplies
- Made in durable black or natural oak finish
- Offered with (4) 3" diameter easy-rolling casters
- Dimensions: 52.5" High, 34" Wide & 19" Deep
- Weights and bands sold separately
- Made in the USA

720631-001	Natural Oak (Weights Sold Separately)	\$699.00
720632-001	Black (Weights Sold Separately).....	699.00
720750-001	Hooks for Weight Carts, 12 Pack	11.95



Tru™ All-In-One Mobile Weight/Storage Rack

- Beautiful craftsmanship on this multi-purpose weight/storage rack
- Almond coated pegboard on front of rack for (9) therapy weights, includes (11) pegboard safety hooks for therapy weights
- Includes (10) dumbbell hooks on left side and (10) hooks on the right side to store dumbbells, can hold (2) 1 to 10 lb. dumbbell sets
- The opposite side of the rack offers three large storage shelves, storage bin for weight bars, miscellaneous bin & resistive band storage area
- Made in durable black or natural oak finish
- Offered with (4) 3" diameter easy-rolling casters
- Dimensions: 66" High, 32" Wide & 24" Deep
- Weights and bands sold separately
- Made in the USA 🇺🇸 🚚

720702-001	Natural Oak (Weights Sold Separately)	\$927.50
720702-002	Black (Weights Sold Separately).....	927.50
720750-001	Hooks for Weight Carts, 12 Pack	11.95



Tru™ Nested Climbing Stools

- Set of 4 stools which nest for storage
- All nested climbing stools are made of ¾" oak plywood
- Wood surfaces are finished with durable clear polyurethane
- Anti-slip treads standard on all platforms to ensure safety
- Each set of stools includes one stool in each of the following dimensions:
- 2"H x 10.5"W x 12"L
- 4"H x 12.5"W x 14"L
- 6"H x 14.5"W x 17.75"L
- 8"H x 16.5"W x 21"L
- Also offered individually
- Weight capacity: 350 lbs.
- Made in the USA

720105-001	Nested Climbing Stools (set of 4).....	\$269.00
720105-002	Nested Climbing Stool 2"	64.95
720105-004	Nested Climbing Stool 4"	69.95
720105-006	Nested Climbing Stool 6"	74.95
720105-008	Nested Climbing Stool 8"	79.95



Tru™ 4-In-1 Weights/Band/Bars Mobile Weight/Storage Rack

- Beautiful craftsmanship on this multi-purpose weight/storage rack
- Almond coated pegboard on front of rack for (9) therapy weights, includes (11) pegboard safety hooks for therapy weights
- Includes (20) dumbbell hooks on the opposite side to store dumbbells, can hold (2) 1 to 10 lb. dumbbell sets
- The left side of the rack offers storage for resistive exercise band, on right side of rack offers storage for weighted exercise bars
- Made in durable black or natural oak finish
- Offered with (4) 3" diameter easy-rolling casters
- Dimensions: 59" High, 30" Wide & 18" Deep
- Weights and bands sold separately
- Made in the USA 🇺🇸 🚚

720747-001	Natural Oak (Weights Sold Separately)	\$735.00
720747-002	Black (Weights Sold Separately).....	735.00
720750-001	Hooks for Weight Carts, 12 Pack	11.95



Tru™ Deluxe Oak Mobile Posture Mirror

- Hand-crafted with high quality solid oak
- All wood surfaces are finished with durable clear polyurethane
- Offers .25" distortion-free mirror
- Offered with (4) easy-rolling casters
- Glass size is 23" x 67", overall size is 28" x 79"
- Made in the USA 🇺🇸 🚚

720169-001	Deluxe Oak Posture Mirror	\$475.00
------------	---------------------------------	----------





Tru™ Deluxe Lift & Place Rack Station

- Hand-crafted with high quality oak plywood
- Wood surfaces are finished with durable clear polyurethane
- Designed to simulate lifting and loading tasks
- Two adjustable oak plywood shelves insert into side mounted rails
- Shelves adjust from 28" to 59" high
- Each adjustable shelf has 100 lb. weight capacity
- Rack station must be secured to wall
- Size is 30"W x 18"D x 78"H
- Made in the USA

720665-001 Lift & Place Rack Station \$805.00



Tru™ Wall-Mount Exercise Ball Rack

- A great storage solution for your roaming exercise balls
- Hand-crafted with high quality solid oak
- All surfaces are finished with durable clear polyurethane
- Available in single exercise ball model or four exercise ball model
- Single exercise ball model dimensions; 5"H x 19.25"W x 24"D
- Four exercise ball model dimensions; 5"H x 66"W x 24"D
- Installation hardware not included
- Made in the USA

720810-001 Multiple Exercise Ball Rack..... \$95.00

720810-002 Single Exercise Ball Rack.....47.50



Tru™ Therapy Putty

- Ideal for increasing hand muscle strength and improving range of motion
- Available in 2 oz., 3 oz., 4 oz., 6 oz., 1 lb. and 5 lb. sizes
- Color coded to reflect firmness: Tan (X-Soft), Yellow (Soft), Red (Med/Soft), Green (Med/Firm) and Blue (Firm)
- Latex Free and non-toxic
- Made in the USA



COLOR CHART

(010)	Tan/X-Soft
(015)	Yellow/Soft
(020)	Red/Med-Soft
(025)	Green/Med-Firm
(030)	Blue/Firm

Therapy Putty

720191-(***)	2 oz		\$3.60
720196-(***)	3 oz		5.05
720192-(***)	4 oz		6.10
720193-(***)	6 oz		8.45
720194-(***)	16 oz		21.50
720195-(***)	5 lb		81.60
720191-(2**)	2 oz Private Label		3.60
720196-(2**)	3 oz Private Label		5.05
720192-(2**)	4 oz Private Label		6.10
720193-(2**)	6 oz Private Label		8.45
720194-(2**)	16 oz Private Label		21.50
720195-(2**)	5 lb Private Label		81.60
720191-100	Private Label Set Up		40.00

Plastic Container Cups

720196-001	Cups 2 & 3 oz		\$1.50
720196-002	Cups 4 & 6 oz		2.25
720196-003	Cups 16 oz		2.50
720196-004	Cups 5 lb		4.25



Cushions & Rolls



Customize Tru™ Brand Products With Your Clinic Name & Logo

A minimum of 36 Tru "Signature Series" Cervical Rolls, Lumbar Rolls or Back Cushions are required. For Silk-screening or private labeling you may mix and match products and sizes to meet the minimum quantity. There is a one-time Setup Charge of \$80.00 (720260-001) with "camera-ready artwork" for the customized name and logo. Logo dimension size is 3" x 5".

Tru™ Lumbar Rolls

The human spine is naturally curved to provide support to the body, however while sitting in traditional chairs, the spine becomes slumped adding stress to the lumbar back. Tru™ Lumbar Rolls increase support to the lumbar spine and assist in maintaining good posture.

- Offered in two popular styles; full size (12"L x 5" Dia.) and half size (12"L x 3"D x 6.25"H)
- Made from high-density urethane foam covered with a durable poplin cover
- Provides standard and firm support and comfort for the lower back
- Promotes proper posture while sitting
- Offered with or without an attachment strap
- Made in the USA

Standard Lumbar Rolls

720227-001	Standard	\$11.75
720227-002	Standard, w/Strap	12.75
720227-003	Standard w/Silk Screen	12.75
720227-004	Standard w/Strap & Silk Screen	13.75
720260-001	Silk Screen Set Up Fee	89.00

Firm Lumbar Rolls

720228-001	Firm	\$11.75
720228-002	Firm w/Strap	12.75
720228-003	Firm w/Silk Screen	12.75
720228-004	Firm w/Strap & Silk Screen	13.75
720260-001	Silk Screen Set Up Fee	89.00



Signature Series

Half Lumbar Rolls

720226-001	Half	\$11.75
720226-002	Half w/Strap	12.75
720226-003	Half w/Silk Screen	12.75
720226-004	Half w/Strap & Silk Screen	13.75
720260-001	Silk Screen Set Up Fee	89.00

Tru™ Back Cushions

The human spine is naturally curved to provide support to the body, however while sitting in traditional chairs, the spine becomes slumped adding stress to the lumbar back. Tru™ Back cushions increase support to the lumbar spine and assist in maintaining good posture.

- Offered in three popular styles; Standard, Bucket Seat and Hugger
- Made from high-density urethane foam covered with a durable poplin cover
- Provides firm support and comfort for the lower back
- Promotes proper posture while sitting
- Offered with or without an attachment strap
- Made in the USA

Signature Series



Standard Back Cushion

720231-001	Cushion	\$20.95
720231-002	w/Strap	21.95
720231-003	Cushion w/Silk Screen	21.95
720231-004	Cushion w/Strap & Silk Screen	22.95
720260-001	Silk Screen Set Up Fee	89.00

Bucket Seat Back Cushion

720233-001	Cushion	\$22.95
720233-002	w/Strap	23.95
720233-003	Cushion w/Silk Screen	23.95
720233-004	Cushion w/Strap & Silk Screen	24.95
720260-001	Silk Screen Set Up Fee	89.00

Hugger Back Cushion

720232-001	Cushion	\$22.95
720232-002	Cushion w/Strap	23.95
720232-003	Cushion w/Silk Screen	23.95
720232-004	Cushion w/Strap & Silk Screen	24.95
720260-001	Silk Screen Set Up Fee	89.00

Tru™ Cervical Roll Insert

Signature Series

- Provides firm support and comfort for the cervical spine
- Made from high-density urethane foam covered with a durable poplin cover
- Offered in standard and firm support
- Designed for use with patient's existing sleeping pillow
- Measures 20" Long x 3.5" Diameter

- Made in the USA
- Check out Tru Signature Series Information on page 7

720219-001	Standard	\$12.50
720220-001	Firm	13.50
720219-003	Standard w/ Silk Screen	13.50
720220-003	Firm w/ Silk Screen	14.50
720260-001	Silk Screen Set Up Fee	89.00



Tru™ Cervical Pillow

- Cervical pillow most recommended by professionals
- Two lobes—one small, one large contour to fit every size neck
- Three foam densities available; Soft, Medium & Firm
- Full 22" length fits most pillowcases
- 18" travel size also available
- Made in the USA

720217-001	22" Soft	24.50
720217-002	22" Medium	24.50
720217-003	22" Firm	24.50



Tru™ Health Solutions Kiosk

- Allows you to display great patient products
- Available in a 36" tall two shelf kiosk and a 72" tall four shelf kiosk
- You can put your company logo on the banner for a low setup charge
- A great way to promote the products you want your customers to use
- Made in the USA

720714-001	2 Shelf Kiosk, 24"W x 36"H x 14"D	\$30.00
720715-001	4 Shelf Kiosk, 24"W x 72"H x 14"D	60.00
720715-100	Private Label Set-Up Fee	50.00

Tru™ Hi-Low Stools

Tru™ Standard Pneumatic Hi-Low Stool

- Hand-operated pneumatic height adjustments from 18.25" to 23.5"
- 5 leg base for increased stability
- Full 16" diameter seat
- 4" medium-density foam cushion for comfort
- Standard: 250 lb. distributed load weight capacity
- Offered in 8 stocking Naugahyde® colors 



 [Color Chart Page 14](#)

720765-002	Standard Pneumatic Hi-Lo Stool, Black	\$85.00
720765-(***)	All Other Colors	92.00

Tru™ Deluxe Pneumatic Hi-Low Stool

- Hand-operated pneumatic height adjustments from 18.25" to 23.5"
- 5 leg base for increased stability
- Full 16" diameter seat
- 4" medium-density foam cushion for comfort
- Deluxe: 300 lb. distributed load weight capacity
- Offered in 8 stocking Naugahyde® colors 



 [Color Chart Page 14](#)

720245-002	Deluxe Pneumatic Hi-Lo Stool, Black	\$123.50
720245-(***)	All Other Colors	130.50



Tru™ Deluxe Training Stairs

- Stairs are hand-crafted with high quality oak plywood
- Includes solid oak handrails and uprights
- All wood surfaces are finished with durable clear polyurethane
- Handrails are ergonomically designed for comfortable grip and ease of use
- Offers (2) 6" high steps and (3) 4" high steps. Steps are 30" wide
- Platform size is 30"W x 24"D x 12"H
- Anti-slip treads standard on all steps and platform to ensure safety
- Weight capacity: 400 lbs.
- Made in the USA  

720515-001	Deluxe Training Stairs	\$945.00
720515-002	Deluxe Training Stairs With Bus Step	1,025.00

Training Stairs

**Tru™ Deluxe Closed-End Staircase**

- Stairs are hand-crafted with high quality oak plywood
- Includes solid oak handrails and uprights
- All wood surfaces are finished with durable clear polyurethane
- Handrails are ergonomically designed for comfortable grip and ease of use
- Handrails enclose platform
- Offers (3) 6" high steps. Steps are 30" wide
- Platform size is 30"W x 24"D x 18"H
- Anti-slip treads standard on all steps and platform to ensure safety
- Weight capacity: 350 lbs.
- Made in the USA  

720639-001 Deluxe Closed-End Staircase \$1,025.00

Tru™ Deluxe 4-Step Closed-End Staircase

- Staircase is hand-crafted with high quality oak plywood
- Includes solid oak handrails and uprights
- All wood surfaces are finished with durable clear polyurethane
- Handrails are ergonomically designed for comfortable grip and ease of use
- Handrails enclose platform
- Offers (4) 6" high steps. Steps are 30" wide
- Platform size is 30"W x 24"D x 24"H
- Requires 9 ft. minimum ceiling
- Anti-slip treads standard on all steps and platform to ensure safety
- Weight capacity: 350 lbs.
- Made in the USA  

720870-001 Tru Deluxe 4-Step Closed-End Staircase \$1,095.00



◀ Corner
Staircase
Setup



▶ Straight
Staircase Setup

Tru™ Deluxe Convertible Staircase

- Stairs are hand-crafted with high quality oak plywood
- Includes solid oak handrails and uprights
- All wood surfaces are finished with durable clear polyurethane
- Designed to setup in straight line or corner configurations
- 11" deep yellow nosed steps enhance safety and durability
- Offers (3) 6" high steps on one end, and (4) 4 ½" high steps on other end
- Platform is 36" x 36" x 18" high
- Anti-slip treads standard on all steps and platform to ensure safety
- Handrails are ergonomically designed for comfortable grip and have 12" ADA overhang
- Weight capacity: 350 lbs.
- Meets ADA guidelines
- Made in the USA  

720556-001 "ADA Compliant" Convertible Staircase \$2,725.00


Tru™ Adj. Height & Width Parallel Bars

- Width adjusts from 15" to 28" with ergonomic control knobs on each upright
- Each upright telescopes up and down in 1 ½" increments and locks into (10) height positions with fail-safe ball-tip locking pin
- Adjustment range is 29" to 42"

- 30" clearance between uprights
- Offered in a satin-finish oak platform with tapered ends for easy wheelchair access and floor mount
- (2) "anti-slip" treads on each end for increased safety
- 1 5/8" diameter one-piece stainless steel handrails
- Heavy gauge white powder coated steel uprights and fittings
- Offered in 7 ft. and 10 ft. lengths
- Platform is 38" wide
- Weight capacity is 500 lbs.
- Made in the USA  

720759-001	7 ft. Floor Mount Parallel Bars.....	\$915.00
720760-001	10 ft. Floor Mount Parallel Bars.....	965.00
720761-001	7 ft. Platform Mounted Parallel Bars.....	1,525.00
720762-001	10 ft. Platform Mounted Parallel Bars.....	1,850.00


Neckpro® II OTD Traction Device

- Simple, easy, effective, portable cervical traction
- No cumbersome weights, water bags or pulleys
- Aids in returning normal neck mobility, flexibility and range of motion
- The ratcheting mechanism enables incremental adjustment of the tension without having to dismantle and reassemble the unit
- Tension is referenced by an audible clicking sound of the ratchet mechanism.
- Maximum Traction pull is 20 lbs.
- Made in the USA

730478-001	NeckPro II OTD Traction Device	\$72.50
730479-001	NeckPro II Head Halter, Each	20.00

Treatment Tables


 Shown With
Optional Lift Back

Table Features Hand-Crafted
High Quality Solid Oak

 Color Chart Page 14

Tru™ Deluxe Oak Treatment Table

- Hand-crafted with high quality solid oak
- All wood surfaces are finished with durable clear polyurethane
- Table top includes 2" medium-density foam and durable Naugahyde® upholstery
- "H" brace is standard on all models for extra stability
- Many options are available to customize your treatment table just the way you want
- Available in 27" W x 72" L, 30" W x 72" L and 30" W x 78" L
- Surface height is 30"
- Offered in 8 stocking Naugahyde® colors
- Weight Capacity: 500 lbs.
- Made in the USA  

730305-(***)	27" X 72"	\$535.00
730001-(***)	30" X 72"	535.00
730002-(***)	30" X 78"	550.00

Tru™ Standard Oak Treatment Table



Features:

Hand-Crafted
Oak Plywood
& Solid Oak Legs

Quality & Simplicity

For multiple treatment applications.

- All wood surfaces are finished with durable clear polyurethane
- Table top includes 1½" medium-density foam and durable Naugahyde® upholstery
- "H" brace is standard on all models for extra stability
- Many options are available to customize your treatment table just the way you want



▶ See Table Accessories Pg. 15

- Available in 27" W x 72" L, 30" W x 72" L and 30" W x 78" L
- Surface height is 30"
- Offered in 8 stocking Naugahyde® colors
- Weight Capacity: 400 lbs.
- Made in the USA

730306-(***) 27" X 72"	\$455.00
730003-(***) 30" X 72"	455.00
730004-(***) 30" X 78"	470.00

Tru™ Heavy-Duty Oak Treatment Table

- Hand-crafted with high quality solid oak
- All wood surfaces are finished with durable clear polyurethane
- Table top includes 2" medium-density foam and durable Naugahyde® upholstery
- Six solid oak legs and a "H" brace is used for extra stability
- Many options are available to customize your treatment table just the way you want
- Available in 30" W x 78" L
- Surface height is 30"
- Offered in 8 stocking Naugahyde® colors
- Weight Capacity: 700 lbs.
- Made in the USA

730354-(***) 30" x 78" \$655.00



▶ Shown With Optional
Lift Back & Divided Shelf



877-882-7844
www.tru-medical.com

Tru™ Upholstery Color Chart
(Additional charge may be incurred for custom colors) To specify color add the appropriate color code after item number (ie. 123456-(***) in Imperial Blue order 123456-003)

002 BLACK	003 IMPERIAL BLUE	004 DOVE GRAY	005 BURGUNDY
008 SAND	012 ROYAL BLUE	052 GUN METAL	053 CLARET
099 CUSTOM	Table Wood Colors		
	NATURAL	MEDIUM BROWN	BLACK


Tru™ Treatment Table Options

We offer 12 different treatment table options to accessorize your Tru™ brand treatment tables.

- 16" Naugahyde® brand. Upholstered shelf is ideal for storage.
- Oak plywood shelf is made from oak plywood to match the wood grain of your table
- Nose hole adds comfort to patients during prone positioning
- Nose hole plug offers continuity to table surface when nose hole is not in use
- Lift-Back option is fully adjustable from 0° to 85°
- Divided shelf provides convenient stackable storage with two compartments
- Paper holder offers convenient storage of table paper
- Paper cutter is constructed of clear plastic, cuts paper evenly
- Soft Touch 3" Dual layer foam offers ultimate comfort
- Edge guard helps protect your tables from damage.
- Made in the USA

Table Options

A	730007-(***) Upholstered Shelf 16" Wide	\$36.95
B	730008-001 Oak Plywood Shelf	55.50
C	730010-001 Nose Hole	33.25
D	730010-002 Nose Hole Plug	28.95
E	730011-001 Lift Back	173.50
F	730012-001 16" Divided Shelf	127.50
G	730013-001 19" x 19" Drawer	107.50
H	730416-001 Modality Drawer, 16"W x 15"D	99.95
I	730537-001 Paper Cutter	17.95
J	730532-001 Paper Holder	39.95
K	730540-001 Soft Touch Foam	75.00
L	730541-001 Edge Guard, Per Foot	3.95



High Quality Accessories
For Your Every Need



Hi-Low Treatment Table



Tru™ Electric Hi-Low Treatment Table

- Hand-crafted with high quality solid oak and 120V, 60 Hz motor that adjusts the height
- Electronically-controlled height adjustment from 22" to 34"
- Comes standard with hand hi-low controller
- All wood surfaces are finished with durable clear polyurethane
- Table top includes 2" medium-density foam and durable Naugahyde® upholstery
- Available in 30" W x 78" L
- Hand control is standard
- Offered in 8 stocking Naugahyde® colors
- Weight Capacity: 500 lbs.
- Made in the USA



730369-(***) 30" x 78"	\$2,189.00
730286-001 Optional Hi-Low Foot Switch	215.00
730286-002 Optional Hi-Low Hand Switch	105.00

Color Chart Page 14

Hand Therapy Table



Tru™ Deluxe Laminate Top Hand Therapy / Work Table

- This Tru™ deluxe work table is great for hand therapy
- Legs & Skirt hand-crafted of solid oak
- Top is made of a smooth laminate material
- Dimensions: 20"L x 42"W x 28"H
- Made in the USA

730525-001 Hand Therapy/Work Table, 20" x 42" \$425.00

Tru™ 3 Section Hi-Low Treatment Table

The Tru™ 3-Section Hi-Low Treatment Table is built to last. The frame is built of high quality square section steel tubing for maximum rigidity and stability. The Tru™ table offers an attractive closed base system that hides all electrical components making the table look spotless. The Tru™ 3-Section Hi-Low Treatment Table exhibits functionality, refinement and great looks.

- Overall Length: 77.75"L
- Overall Width: 26"W
- Head Section: 16.25"L x 26.5"W
- Middle Section: 26"L x 26.5"W
- Foot Section: 34.75"L x 26.5"W
- Max Height: 37"
- Min Height: 17.75"
- Hi-Low Controls: Foot or All-Round Bar Control (ARB)
- 2" Firm Density Foam Tops
- Offered in (8) Upholstery Colors
- Integrated, Retractable 3" Casters
- 110 Volts, 50-60Hz Motor
- 500 lb. Lifting Weight Capacity
- Table Weight: 210 lbs.



730459-(***) Table w/Foot Switch, Black & Blue \$1,895.00
 730459-(***) Table w/Foot Switch, All Other Colors 1,995.00
 730460-(***) Table w/Foot Bar, Black & Blue 2,260.00
 730460-(***) Table w/Foot Bar, All Other Colors 2,360.00



Color Chart Page 14



Shown With
Optional Shelf

Tru™ Trainers Taping Table

- Hand-crafted with high quality solid oak
- All wood surfaces are finished with durable clear polyurethane
- Table top includes 2" medium-density foam and durable Naugahyde® upholstery
- Available in 24"W x 48"L x 36"H
- Can also choose an optional shelving system for storage
- Offered in 8 stocking Naugahyde® colors
- Weight Capacity: 450 lbs.
- Made in the USA

730183-(***) Taping Table..... \$495.00

730184-001 Optional Shelf 115.00



Tru™ Deluxe Oak 2-Seat Taping Station

- A great taping station solution for the most demanding training rooms
- (2) Naugahyde® upholstered removable taping modules are 36"H x 42"D x 26"W
- These sections offer a front-facing storage shelf, a 2-shelf outside facing storage solution and a storage cabinet with doors that have concealed hinges
- These modules are separated by a center section constructed of ¾" oak plywood that is 32"high
- This section offers a drawer and storage cabinet with doors that have concealed hinges
- Offered in 8 stocking Naugahyde® colors
- Made in durable black or natural oak finish
- Accessories offered are a trash can, drawer lock and door lock
- Custom configurations of this station can be offered per request
- Made in the USA

730397-(***) 6ft. L x 42" D x 36" H..... \$2,405.00

730398-001 Hamper/Trash-Can Option for Taping Station 130.00

730462-001 Optional Lock for Double Cabinet Door..... 79.95

730463-001 Optional Top Drawer Lock 45.00

Mat Tables



► Color Chart Page 14

Tru™ Deluxe Oak Mat Table

- Hand-crafted with high quality solid oak
- All wood surfaces are finished with durable clear polyurethane
- Table top includes 2" medium-density foam and durable Naugahyde® upholstery
- Offered in three sizes, 4' x 7', 5' x 7' and 6' x 8'
- Surface height is 20"
- Offered in 8 stocking Naugahyde® colors
- Weight Capacity: 500 lbs.
- Made in the USA 

730057-(***) 4' x 7'	\$710.00
730058-(***) 5' x 7'	940.00
730059-(***) 6' x 8'	1085.00



► Color Chart Page 14

Tru™ Standard Oak Mat Table

- Hand-crafted with high quality oak plywood & solid oak legs
- All wood surfaces are finished with durable clear polyurethane
- Table top includes 1½" medium density foam and durable Naugahyde® upholstery
- Offered in three sizes, 4' x 7', 5' x 7' and 6' x 8'
- Surface height is 20"
- Offered in 8 stocking Naugahyde® colors
- Weight Capacity: 400 lbs.
- Made in the USA 

730060-(***) 4' x 7'	\$620.00
730061-(***) 5' x 7'	830.00
730062-(***) 6' x 8'	985.00



*Accessories Not Included

Tru™ "Low Rider" Storage Tray For Mat Tables

This mobile storage tray is a great idea for the clinic with little to no storage area

- Simple put storable items in the tray and store under your mat table
- Hand-crafted with high quality oak plywood
- All wood surfaces are finished with durable clear polyurethane
- Includes 1" diameter easy-rolling casters
- Offered in two sizes; 21" x 24" and 36" x 42"
- Overall height is 7"
- Weight Capacity: 100 lbs.
- Made in the USA

730395-001 Storage Tray, 21" x 24"	\$190.00
730395-002 Storage Tray, 36" x 42"	215.00

Work Table



Tru™ Deluxe Hi-Low Butcher Block Work Table

- Great solution for work/hand therapy
- Offered in two Hi-Low options; electric and manual crank
- Top hand-crafted of quality 1.5" solid hardwood
- Top is finished with durable clear polyurethane
- Top dimensions are 36" x 60"
- Height adjusts from 27" to 38"
- Made in the USA 

730400-001 Electric Hi-Low Table, 36" x 60"	\$3,095.00
730400-002 Crank Hi-Low Table, 36" x 60"	2,675.00



Down Position



Features:

Hand-crafted high quality solid oak

Space Saving Folded Position

Color Chart Page 14



Color Chart Page 14

Tru™ Wall-Mounted Folding Mat Table

- All wood surfaces are finished with durable clear polyurethane
- Table top includes 2" medium-density foam and durable Naugahyde® upholstery
- Offered in two sizes, 3' x 7' and 4' x 7'
- Folds up to the wall when not in use, saving valuable floor space
- Surface height in unfolded position is 20"
- Offered in 8 stocking Naugahyde® colors
- Weight Capacity: 500 lbs.
- Made in the USA 

730366-(***) 3' x 7' Folding Mat Table	\$980.00
730066-(***) 4' x 7' Folding Mat Table	980.00

Tru™ Deluxe Split Lift Back Mat Table

- Hand-crafted with high quality solid oak
- All wood surfaces are finished with durable clear polyurethane
- Table top includes 2" medium-density foam and durable Naugahyde® upholstery
- Offers unique split lift-back sections to help position the patient with infinite adjustment from 0 degrees to 85 degrees
- Offered in two sizes, 4' x 7' and 5' x 7'
- Surface height is 20"
- Offered in 8 stocking Naugahyde® colors
- Weight Capacity: 500 lbs.
- Made in the USA 

730204-(***) 4' x 7' Table	\$1,210.00
730244-(***) 5' x 7' Table	1,415.00

Electric Hi-Low Mat Table



Color Chart Page 14

OFFERED IN THREE SIZES AND 8 COLORS

730526-(***) Manual Crank 4' x 7'	\$2,795.00
730527-(***) Manual Crank 5' x 7'	2,895.00
730528-(***) Manual Crank 6' x 8'	3,150.00

Tru™ Electric Hi-Low Mat Table

- Hand-crafted with high quality solid oak and 120V, 60 Hz motor that adjusts the height
- Electronically-controlled hand switch adjusts table from 22" to 34"
- Table top includes 2" medium-density foam and durable Naugahyde® upholstery
- Offered in three sizes, 4' x 7', 5' x 7' and 6' x 8'
- Hand control is standard
- Offered in 8 stocking Naugahyde® colors
- Weight capacity 700 lbs.
- Made in the USA 



730258-(***) 4' x 7' Electric Hi-Low Mat Table	\$3,089.00
730259-(***) 5' x 7' Electric Hi-Low Mat Table	3,189.00
730260-(***) 6' X 8' Electric Hi-Low Mat Table	3,465.00
730286-001 Optional Foot Switch	215.00
730286-002 Optional Hand Switch	105.00



▶ Color Chart Page 14

Tru™ Electric Hi-Low Split Lift-Back Mat Table

- Hand-crafted with high quality solid oak and 120V, 60 Hz motor that adjusts the height
- Electronically-controlled hand switch adjusts table height from 22" to 34"
- Table top includes 2" medium-density foam and durable Naugahyde® upholstery
- Offers unique split lift-back sections to help position the patient with infinite adjustment from 0 degrees to 85 degrees
- Offered in two sizes, 4' x 7' and 5' x 7'
- Hand control is standard
- Offered in 8 stocking Naugahyde® colors
- Weight capacity 500 lbs.
- Made in the USA 

730350-(***) 4' x 7' Table	\$3,450.00
730351-(***) 5' x 7' Table	3,595.00
730286-001 Optional Foot Switch	215.00
730286-002 Optional Hand Switch	105.00


Bolsters & Wedges



 **Tru™**



 **Tru™**



Tru™ Bolsters

- High quality bolsters offer great versatility in positioning patients
- Improve weight-bearing, dynamic gross motor movements, coordination, righting reactions, balance, strength and muscle tone
- Made with high density urethane foam and covered with durable Naugahyde® upholstery
- Offered in 8 stocking Naugahyde® colors
- Made in the USA

730140-(***) 5" x 12"	\$37.75
730141-(***) Half Bolster, 8" x 16"	59.95
730218-(***) 4" x 24"	41.25
730143-(***) 6" x 24"	51.50
730142-(***) 8" x 18"	49.95
730144-(***) 8" x 24"	65.95
730145-(***) 8" x 36"	96.95
730146-(***) 10" x 36"	116.95
730441-(***) 12" x 36"	157.95
730442-(***) 12" x 48"	198.50
730443-(***) 14" x 48"	229.50
730444-(***) 16" x 48"	270.50

Tru™ Speciality Half Bolsters

- High quality specialty half bolsters offer great versatility in knee flexion and extension
- Improve weight-bearing, dynamic gross motor movements, coordination, righting reactions, balance, strength and muscle tone
- Made with high density urethane foam and covered with durable Naugahyde® upholstery
- Available in a variety of sizes
- Offered in 8 stocking Naugahyde® colors
- Made in the USA

730446-(***) 2" x 4" x 18"	\$48.95
730447-(***) 3" x 4" x 18"	50.95
730448-(***) 4.5" x 6" x 24"	69.95
730449-(***) 6" x 8" x 30"	95.00


Tru™ Wedges

- High quality wedges offer great versatility in positioning patients
- Flat surfaced positioning shapes used to build head control, sitting balance, trunk control and upper extremity function
- Made with high density urethane foam and covered with durable Naugahyde® upholstery
- Available in a variety of sizes
- Offered in 8 stocking Naugahyde® colors
- Made in the USA

730150-(***) 4" x 20" x 22" Wedge.....	\$69.50
730151-(***) 6" x 20" x 22" Wedge.....	82.50
730445-(***) 6" x 24" x 28" Wedge.....	113.50
730152-(***) 8" x 20" x 22" Wedge.....	96.95
730153-(***) 8" x 24" x 28" Wedge.....	113.95
730154-(***) 10" x 20" x 22" Wedge.....	104.00
730155-(***) 10" x 24" x 28" Wedge.....	122.50
730156-(***) 12" x 24" x 28" Wedge.....	135.50



► Color Chart Page 14



Tru™ Specialty Wedges

- High quality specialty wedges offer great versatility in positioning patients
- Made with high density urethane foam and covered with durable Naugahyde® upholstery
- Three styles available; knee crescent, knee wedge and flat-top wedge
- Offered in 8 stocking Naugahyde® colors
- Made in the USA

Tru™ Flat Top Wedge

- Dimensions: 20" W x 12" H x 32" L

730149-(***) Flat Top Wedge\$157.50

Tru™ Knee Wedge

- Dimensions: 20" W x 10½" H x 22" L

730148-(***) Knee Wedge.....\$113.00

Tru™ Knee Crescent

- Dimensions: 13" W x 9" H x 20" L


730147-(***) Knee Crescent \$96.00

Tru™ Anti-Microbial Positioning Bed Wedge

- Great solution for patients who need proper support while sleeping
- Offered in a 30° and 45° wedge angles
- Made with high density urethane foam
- Covered in a black anti-Microbial fabric that is durable and soft to the touch
- Made in the USA

730483-001 7" x 12" x 24" Black \$73.95

730484-001 8" x 10" x 28" Black..... 75.95



► Color Chart Page 14



► Color Chart Page 14

Tru™ Naugahyde® Pillows

- Great solution for the busy clinic
- Offered in two different sizes of 12" x 15" and 15" x 30"
- No pillowcase needed, just disinfect between patients
- Filled with hypo-allergenic polyester fill
- Offered in 8 stocking Naugahyde® colors
- Made in the USA

730412-(***) 12" x 15" Naugahyde Pillow \$39.95

730413-(***) 15" x 30" Naugahyde Pillow 49.95

Tru™ Naugahyde® Prone Face Pillow

- A great low-cost solution for treating patients in the prone position
- Filled with hypo-allergenic polyester fill
- Offered in 8 stocking Naugahyde® colors
- Dimensions are 15" x 17"
- Made in the USA

730489-(***) Prone Face Pillow 15" x 17" \$30.95

Hook & Loop Fasteners



Measures Body Composition Using Electrodes.

TruBIA™ Body Composition Analyzer

A great new product that measures impedance for body composition using surface electrodes. Offers accurate and repeatable results that will articulate total water analyzes, hydration and wasting analysis. Competitively priced.

- Intercellular and extracellular measurements
- Whole and segmental body composition capability
- Phase angle measurements
- Touch screen interface
- Rechargeable battery
- Stores patient historical data
- Capable of generating reports

770177-001	TruBIA™ Body Composition Analyzer.....	\$2,195.00
770177-002	Electrodes, 100/Package.....	19.95
770177-003	Lead-Wire, Pair	70.00

TruBIA™ Features

• Capacity (pounds)	No Scale – Up to 600 lbs. input
• Position of patient for measurement	Supine
• Contact Type	Surface Electrodes
• Multi-Frequency	Yes
• Full Body Measurements	Yes
• Reactance/Resistance/Phase Angle	Yes
• Segmental Measurement	Yes
• History	Yes
• External Print Capability	Yes
• USB Connection	Yes
• DEXA Comparison	No
• Air Displacement Comparison	+/- 2%

Tru™ Hook & Loop Fasteners

- A great solution for hand therapy splinting creations
- These fasteners consist of two mating components: hook & loop
- Choose between sew-on hook and loop or pressure sensitive hook and loop
- Made of durable nylon material
- Offered in three different colors: White, Beige & Black
- Offered in 1" width and 2" width
- Offered in 25 yard rolls or 10 yard rolls



Color Chart


(001)-White



(002)-Black



(003)-Beige



Hook

750651-(***)	1" Hook 10 yd.....	\$9.95
750652-(***)	1" Hook 10 yd. Pressure Sensitive	16.50
750653-(***)	1" Hook 25 yd.....	21.95
750654-(***)	1" Hook 25 yd. Pressure Sensitive	36.95
750655-(***)	2" Hook 10 yd.....	17.95
750656-(***)	2" Hook 10 yd. Pressure Sensitive	29.95
750657-(***)	2" Hook 25 yd	36.95
750658-(***)	2" Hook 25 yd. Pressure Sensitive.....	59.95

Loop

750659-(***)	1" Loop 10 yd	\$9.95
750660-(***)	1" Loop 10 yd. Pressure Sensitive.....	16.50
750661-(***)	1" Loop 25 yd	21.95
750662-(***)	1" Loop 25 yd. Pressure Sensitive	36.95
750663-(***)	2" Loop 10 yd	17.95
750664-(***)	2" Loop 10 yd. Pressure Sensitive	29.95
750665-(***)	2" Loop 25 yd	36.95
750666-(***)	2" Loop 25 yd. Pressure Sensitive	59.95



Tru™ Elasticated Tubular Bandages

An elasticated tubular bandage used in the treatment of a variety of conditions where tissue support is required.

- Knitted from 100% cotton
- Covered elastic threads are laid into the cotton material to form a continuous spiral, no pins or tape needed
- Offered in 6 sizes for every body part
- Each box contains one 33 ft. roll
- Made in the USA

LATEX FREE

790079-001	Size A, 2½"	\$19.25
790079-002	Size B, 2¾"	28.50
790079-003	Size C, 3"	31.50
790079-004	Size D, 3½"	36.95
790079-005	Size E, 4"	38.95
790079-006	Size F, 5"	41.25
790079-007	G, 5" x 33 ft.	47.95
790079-008	J, 6.75" x 33 ft.	60.95



Tru™ Splint Liners

- Two-ply cotton stockinette (uni-size) liners are thick & soft making for a great hand/wrist splint liner
- Offers finished seams and smooth edges
- Re-usable, machine wash and dry
- Offered in 3" x 6" and 3" x 13" with and without thumb
- Package of 10
- Made in the USA

LATEX FREE

750285-001	3" x 6" w/Thumb	\$27.50
750285-002	3" x 6" w/o Thumb	26.50
750286-001	3" x 13" w/Thumb	31.50
750286-002	3" x 13" w/o Thumb	30.75

Tru™ Stockinette



LATEX FREE

Tru™ Polyester Stockinette

- Made of durable 100% polyester for support that will not shrink when washed
- Excellent undergarment for custom splinting with light compression
- Offered in 1", 2", 3", 4" and 6" widths
- Each box contains one 25 yd. roll
- Made in the USA

790078-001	1" Stockinette, 25 yd box	\$10.75
790078-002	2" Stockinette, 25 yd box	13.50
790078-003	3" Stockinette, 25 yd box	20.50
790078-004	4" Stockinette, 25 yd box	26.25
790078-006	6" Stockinette, 25 yd box	39.50



LATEX FREE

Tru™ Cotton Stockinette


- High-quality 100% cotton stockinette
- Knit tubular bandage stretches up to three times its original size for conformability
- Excellent undergarment for custom splinting and cast fabrication
- Offered in 2", 3", 4" and 6" widths
- Each box contains one 25 yd. roll
- Made in the USA

790122-002	2" Stockinette, 25 yd box	\$13.50
790122-003	3" Stockinette, 25 yd box	20.25
790122-004	4" Stockinette, 25 yd box	25.95
790122-006	6" Stockinette, 25 yd box	38.95

Lumbar Supports



► **Compression Pad**



Tru™ Deluxe Lumbo Sacral Support

- Features 10" back with pocket and compression pad
- Hook and loop closure for ease of application and adjustment
- Offered in sized and universal fit belts
- Offered with four flexible stays for increased support
- Sizes based on hip or waist measurement, whichever is larger
- Made in the USA

4-Stay

750133-200	Small (26"-30")		\$29.95
750133-300	Medium (30"-36")		29.95
750133-400	Large (34"-40")		29.95
750133-500	X-Large (38"-46")		29.95
750133-600	XX-Large (48"-58")		32.50
750133-700	XXX-Large		35.95
750133-900	Universal (26"-38")		29.95



◀ **Front View**



▶ **Back View**

Tru™ Universal Lumbar Support

- Tricot laminate sides with elastic back offer wide range of fittings
- Available in 8" and 10" backs for proper lumbar support
- Hook and loop closure for ease of application and adjustment
- 4 flexible stays for additional support
- Elastic side pull for increased support
- Sizes based on hip or waist measurement, whichever is larger
- Made in the USA

750137-100	8" X-Small		\$23.50
750137-500	8" X-Large		23.50
750137-900	8" Universal		23.50
750138-100	10" X-Small		25.50
750138-500	10" X-Large		25.50
750138-900	10" Universal		25.50



ORDER TODAY
Visit Our Website
www.tru-medical.com



Tru™ Universal Cervical Collar

- Contour-cut for comfort under chin
- Made with medium density urethane foam for proper support
- Hook and loop closure for ease of application and adjustment
- Includes washable stockinette cover
- Made in the USA



750151-250	2½" Height		\$6.75
750151-300	3" Height		6.75
750151-350	3½" Height		6.75

Tru™ Crisscross Lumbar Support




MADE IN THE USA

- Features 9" all-knit elastic design for comfort
- Hook and loop closure for ease of application and adjustment
- Features 4" crisscross elastic for rigid support
- Sizes based on hip or waist measurement, whichever is larger

750136-200	Small		\$19.50
750136-300	Medium		19.50
750136-400	Large		19.50
750136-500	X-Large		19.50
750136-600	XX-Large		20.50
750136-900	Universal		19.50



▲ Velcro Straps



▲ Offered In 3 Styles

Tru™ Black Lumbar Supports

- Offered in sized and universal fit belts
- Belts offered in 6", 8" or 10" back for custom fit
- Hook and loop closure for ease of application and adjustment
- Four flexible stays for additional support
- Sizes based on hip or waist measurement, whichever is larger
- Offered with or without suspenders
- Made in the USA

6" Black Lumbar Support

750399-200 X-Small/Small	\$29.50
750399-500 X-Large	29.50
750399-600 XX-Large	31.50
750399-900 Universal	29.50
750400-200 X-Small/Small, w/Suspenders.....	32.50
750400-500 X-Large, w/Suspenders	32.50
750400-600 XX-Large, w/Suspenders	34.50
750400-900 Universal w/Suspenders.....	32.50

8" Black Lumbar Support

750401-200 X-Small/Small	\$31.50
750401-500 X-Large	31.50
750401-600 XX-Large	33.50
750401-900 Universal	31.50
750402-200 X-Small/Small, w/Suspenders	34.50
750402-500 X-Large, w/Suspenders	34.50
750402-600 XX-Large, w/Suspenders	36.50
750402-900 Universal, w/Suspenders	34.50

10" Black Lumbar Support

750289-200 X-Small/Small	\$33.50
750289-500 X-Large	33.50
750289-600 XX-Large	35.50
750289-900 Universal	33.50
750290-200 X-Small/Small, w/Suspenders.....	36.50
750290-500 X-Large, w/Suspenders	36.50
750290-600 XX-Large, w/Suspenders	38.50
750290-900 Universal, w/Suspenders	36.50

Hand Dynamometers



Dual Scale
Readout
Displays
In Lb. & Kg.

Tru™ Hand Dynamometer

- Designed for reliable grip-force evaluations from 0 to 200 lbs. (0 to 90 kg)
 - Dual scale readout displays in lb. and kg.
 - Peak-hold needle automatically retains the highest reading until reset
 - Handle easily adjusts to one of five grip positions from 1 3/8" to 3 3/8", in half inch increments to accommodate any hand size
 - Includes carrying case
 - Two year warranty
- 770001-001 200lb. Hand Dynamometer \$255.00

Tru™ Hand Evaluation Sets

- 3 Piece Set includes; Hydraulic Hand Dynamometer (0 lbs. - 200 lbs.), Hydraulic Pinch Gauge (0 lbs. - 50 lbs.), and 5 1/2" Stainless Steel Finger Goniometer
 - 7 Piece Set also includes; 3 Point Discriminator, Pinwheel, Finger Motion & Circumference Gauges
 - Peak-hold needle automatically retains the highest reading until reset
 - Includes carry case
 - Two year warranty
- 770003-007 7-Piece Set.....\$615.00
770003-003 3-Piece Set..... 449.00



INCLUDES CARRY CASE

Tru™ 50lb. Pinch Gauge

- Highly accurate measurements of 0 - 45 lbs. or (0 to 20 kg)
- Dual scale readout displays in lb. and kg.
- Peak-hold needle automatically retains the highest reading until reset
- Makes for dependable, repeatable test results for all pinch tests (tip, key, and palmar)
- Includes carry case
- Two year warranty

770002-001 Pinch Gauge\$195.00



Tru™ Cloth Tape Measure

- 5 ft. total length retractable measuring tape
- Measures in inches and centimeters
- Made from Woven cloth tape for durability
- One finger retractable push button control

770080-001 Cloth Tape.....\$2.95

Cushions & Wedges




Tru™ Leg Spacer

- Designed to help relieve pressure and strain on the lower back, hips, knees and ankles while you sleep
- Made from high-density urethane foam covered with a durable poplin cover

- While lying on your side, place the Tru™ Leg Spacer between legs just above the knee
- Helps relax muscles
- Made in the USA

780057-002 10" x 8" x 6"\$13.95



Tru™ Posture Wedge

This is an economical way to help improve posture and ease lower back pain while sitting. The incline helps restore and maintain correct ergonomic seating

- Made from high-density urethane solid foam
- Dimensions: 15½" W x 13½" D, 2½" high at rear, ½" at front
- Helps ease lower back pain while sitting
- Helps restore and maintain correct ergonomic seating
- Made in the USA

780060-001 Posture Wedge\$10.25



TruBIA™ Full Body Cloth Measuring Tape

TruBIA™ Full Body Measuring Tape is a great product to take care of all your circumferential body measurements easily and accurately without assistance!

- Retractable tape measure with pin lock can be used as a circumference tape measure by locking into place and retracts easily with the push of a button.

- Measures up to 60 inches in both standard and metric systems
- Soft, Flexible, and durable vinyl tape with two side units
- Measure to 1/16 inch and/or to 1 millimeter
- One finger retractable push button control

770181-001 Full Body Cloth Measuring Tape, 5 ft.\$6.95

TMS Part #.....Page #	TMS Part #.....Page #	TMS Part #.....Page #	TMS Part #.....Page #	TMS Part #.....Page #
710101-001.....2	710281-001.....4	720191-020.....9	720233-004.....10	730478-001.....13
710101-002.....2	710281-002.....4	720191-025.....9	720260-001.....10	730479-001.....13
710101-003.....2	710281-003.....4	720191-030.....9	720217-001.....11	730306-002.....14
710101-011.....2	710281-004.....4	720196-010.....9	720217-002.....11	730306-003.....14
710101-022.....2	710582-001.....4	720196-015.....9	720217-003.....11	730306-004.....14
710101-033.....2	710582-002.....4	720196-020.....9	720245-002.....11	730306-005.....14
710102-001.....2	710582-003.....4	720196-025.....9	720245-003.....11	730306-008.....14
710102-002.....2	710582-004.....4	720196-030.....9	720245-004.....11	730306-012.....14
710102-003.....2	720005-046.....5	720192-010.....9	720245-005.....11	730306-052.....14
710102-011.....2	720005-047.....5	720192-015.....9	720245-008.....11	730306-053.....14
710102-022.....2	720005-048.....5	720192-020.....9	720245-012.....11	730306-099.....14
710102-033.....2	720005-049.....5	720192-025.....9	720245-052.....11	730003-002.....14
710103-001.....2	720005-051.....5	720192-030.....9	720245-053.....11	730003-003.....14
710103-002.....2	720417-046.....5	720193-010.....9	720245-099.....11	730003-004.....14
710103-003.....2	720417-047.....5	720193-015.....9	720515-001.....11	730003-005.....14
710103-004.....2	720417-048.....5	720193-020.....9	720515-002.....11	730003-008.....14
710103-005.....2	720417-049.....5	720193-025.....9	720714-001.....11	730003-012.....14
710103-006.....2	720417-051.....5	720193-030.....9	720715-001.....11	730003-052.....14
710103-008.....2	720485-046.....5	720194-010.....9	720715-100.....11	730003-053.....14
710103-011.....2	720485-047.....5	720194-015.....9	720765-002.....11	730003-099.....14
710103-022.....2	720485-048.....5	720194-020.....9	720765-003.....11	730004-002.....14
710103-033.....2	720485-049.....5	720194-025.....9	720765-004.....11	730004-003.....14
710103-044.....2	720485-051.....5	720194-030.....9	720765-005.....11	730004-004.....14
710103-055.....2	720486-001.....5	720195-010.....9	720765-008.....11	730004-005.....14
710103-066.....2	720486-002.....5	720195-015.....9	720765-012.....11	730004-008.....14
710103-088.....2	720795-001.....5	720195-020.....9	720765-052.....11	730004-012.....14
710103-100.....2	720050-001.....6	720195-025.....9	720765-053.....11	730004-052.....14
770177-001.....2	720050-002.....6	720195-030.....9	720765-099.....11	730004-053.....14
770177-002.....2	720051-001.....6	720196-001.....9	720556-001.....12	730004-099.....14
770177-003.....2	720051-002.....6	720196-002.....9	720639-001.....12	730354-002.....14
710126-001.....3	720640-001.....6	720196-003.....9	720639-002.....12	730354-003.....14
710126-002.....3	720640-002.....6	720196-004.....9	720759-001.....13	730354-004.....14
710126-003.....3	720711-046.....6	720665-001.....9	720760-001.....13	730354-005.....14
710126-004.....3	720711-047.....6	720810-001.....9	720761-001.....13	730354-008.....14
710126-005.....3	720711-048.....6	720810-002.....9	720762-001.....13	730354-012.....14
710126-011.....3	720711-049.....6	720219-001.....10	730305-002.....13	730354-052.....14
710126-022.....3	720711-051.....6	720219-003.....10	730305-003.....13	730354-053.....14
710126-033.....3	720719-046.....6	720220-001.....10	730305-004.....13	730354-099.....14
710126-044.....3	720719-047.....6	720220-003.....10	730305-005.....13	730007-002.....15
710126-055.....3	720719-048.....6	720226-001.....10	730305-008.....13	730007-003.....15
710126-100.....3	720719-049.....6	720226-002.....10	730305-012.....13	730007-004.....15
710471-001.....3	720719-051.....6	720226-003.....10	730305-052.....13	730007-005.....15
710471-002.....3	720806-001.....6	720226-004.....10	730305-053.....13	730007-008.....15
710471-004.....3	720586-001.....7	720227-001.....10	730305-099.....13	730007-012.....15
710471-005.....3	720588-001.....7	720227-002.....10	730001-002.....13	730007-052.....15
710471-011.....3	720587-001.....7	720227-003.....10	730001-003.....13	730007-053.....15
710471-055.....3	720589-001.....7	720227-004.....10	730001-004.....13	730007-099.....15
710472-001.....3	720631-001.....7	720228-001.....10	730001-005.....13	730008-001.....15
710472-002.....3	720632-001.....7	720228-002.....10	730001-008.....13	730010-001.....15
710472-004.....3	720865-001.....7	720228-003.....10	730001-012.....13	730010-002.....15
710472-011.....3	720105-001.....8	720228-004.....10	730001-052.....13	730011-001.....15
710624-001.....3	720105-002.....8	720231-001.....10	730001-053.....13	730012-001.....15
710624-002.....3	720105-004.....8	720231-002.....10	730001-099.....13	730013-001.....15
710624-003.....3	720105-006.....8	720231-003.....10	730002-002.....13	730532-001.....15
710624-004.....3	720105-008.....8	720231-004.....10	730002-003.....13	730537-001.....15
710065-001.....4	720169-001.....8	720232-001.....10	730002-004.....13	730416-001.....15
710066-001.....4	720702-001.....8	720232-002.....10	730002-005.....13	730536-001.....15
710072-001.....4	720702-002.....8	720232-003.....10	730002-008.....13	730536-100.....15
710280-001.....4	720747-001.....8	720232-004.....10	730002-012.....13	730369-002.....15
710280-002.....4	720747-002.....8	720233-001.....10	730002-052.....13	730369-003.....15
710280-003.....4	720191-010.....9	720233-002.....10	730002-053.....13	730369-004.....15
710280-004.....4	720191-015.....9	720233-003.....10	730002-099.....13	730369-005.....15

TMS Part #.....Page #	TMS Part #.....Page #	TMS Part #.....Page #	TMS Part #.....Page #	TMS Part #.....Page #
730369-008.....15	730058-012.....17	730204-002.....18	730527-099.....18	730144-053.....19
730369-012.....15	730058-052.....17	730204-003.....18	730528-002.....18	730144-099.....19
730369-052.....15	730058-053.....17	730204-004.....18	730528-003.....18	730145-002.....19
730369-053.....15	730058-099.....17	730204-005.....18	730528-004.....18	730145-003.....19
730369-099.....15	730059-002.....17	730204-008.....18	730528-005.....18	730145-004.....19
730183-002.....16	730059-003.....17	730204-012.....18	730528-008.....18	730145-005.....19
730183-003.....16	730059-004.....17	730204-052.....18	730528-012.....18	730145-008.....19
730183-004.....16	730059-005.....17	730204-053.....18	730528-052.....18	730145-012.....19
730183-005.....16	730059-008.....17	730204-099.....18	730528-053.....18	730145-052.....19
730183-008.....16	730059-012.....17	730244-002.....18	730528-099.....18	730145-053.....19
730183-012.....16	730059-052.....17	730244-003.....18	730140-002.....19	730145-099.....19
730183-052.....16	730059-053.....17	730244-004.....18	730140-003.....19	730146-002.....19
730183-053.....16	730059-099.....17	730244-005.....18	730140-004.....19	730146-003.....19
730183-099.....16	730060-002.....17	730244-008.....18	730140-005.....19	730146-004.....19
730184-001.....16	730060-003.....17	730244-012.....18	730140-008.....19	730146-005.....19
730397-002.....16	730060-004.....17	730244-052.....18	730140-012.....19	730146-008.....19
730397-003.....16	730060-005.....17	730244-053.....18	730140-052.....19	730146-012.....19
730397-004.....16	730060-008.....17	730244-099.....18	730140-053.....19	730146-052.....19
730397-005.....16	730060-012.....17	730258-002.....18	730140-099.....19	730146-053.....19
730397-008.....16	730060-052.....17	730258-003.....18	730218-002.....19	730146-099.....19
730397-012.....16	730060-053.....17	730258-004.....18	730218-003.....19	730441-002.....19
730397-052.....16	730060-099.....17	730258-005.....18	730218-004.....19	730441-003.....19
730397-053.....16	730061-002.....17	730258-008.....18	730218-005.....19	730441-004.....19
730397-099.....16	730061-003.....17	730258-012.....18	730218-008.....19	730441-005.....19
730398-001.....16	730061-004.....17	730258-052.....18	730218-012.....19	730441-008.....19
730462-001.....16	730061-005.....17	730258-053.....18	730218-052.....19	730441-012.....19
730463-001.....16	730061-008.....17	730258-099.....18	730218-053.....19	730441-052.....19
730459-002.....16	730061-012.....17	730259-002.....18	730218-099.....19	730441-053.....19
730459-003.....16	730061-052.....17	730259-003.....18	730141-002.....19	730441-099.....19
730459-004.....16	730061-053.....17	730259-004.....18	730141-003.....19	730442-002.....19
730459-005.....16	730061-099.....17	730259-005.....18	730141-004.....19	730442-003.....19
730459-008.....16	730062-002.....17	730259-008.....18	730141-005.....19	730442-004.....19
730459-011.....16	730062-003.....17	730259-012.....18	730141-008.....19	730442-005.....19
730459-012.....16	730062-004.....17	730259-052.....18	730141-012.....19	730442-008.....19
730459-052.....16	730062-005.....17	730259-053.....18	730141-052.....19	730442-012.....19
730459-053.....16	730062-008.....17	730259-099.....18	730141-053.....19	730442-052.....19
730459-099.....16	730062-012.....17	730260-002.....18	730141-099.....19	730442-053.....19
730460-002.....16	730062-052.....17	730260-003.....18	730142-002.....19	730442-099.....19
730460-003.....16	730062-053.....17	730260-004.....18	730142-003.....19	730443-002.....19
730460-004.....16	730062-099.....17	730260-005.....18	730142-004.....19	730443-003.....19
730460-005.....16	730395-001.....17	730260-008.....18	730142-005.....19	730443-004.....19
730460-008.....16	730395-002.....17	730260-012.....18	730142-008.....19	730443-005.....19
730460-011.....16	730400-001.....17	730260-052.....18	730142-012.....19	730443-008.....19
730460-012.....16	730400-002.....17	730260-053.....18	730142-052.....19	730443-012.....19
730460-052.....16	730366-002.....18	730260-099.....18	730142-053.....19	730443-052.....19
730460-053.....16	730366-003.....18	730526-002.....18	730142-099.....19	730443-053.....19
730460-099.....16	730366-004.....18	730526-003.....18	730143-002.....19	730443-099.....19
730525-001.....16	730366-005.....18	730526-004.....18	730143-003.....19	730444-002.....19
730057-002.....17	730366-008.....18	730526-005.....18	730143-004.....19	730444-003.....19
730057-003.....17	730366-012.....18	730526-008.....18	730143-005.....19	730444-004.....19
730057-004.....17	730366-052.....18	730526-012.....18	730143-008.....19	730444-005.....19
730057-005.....17	730366-053.....18	730526-052.....18	730143-012.....19	730444-008.....19
730057-008.....17	730366-099.....18	730526-053.....18	730143-052.....19	730444-012.....19
730057-012.....17	730066-002.....18	730526-099.....18	730143-053.....19	730444-052.....19
730057-052.....17	730066-003.....18	730527-002.....18	730143-099.....19	730444-053.....19
730057-053.....17	730066-004.....18	730527-003.....18	730144-002.....19	730444-099.....19
730057-099.....17	730066-005.....18	730527-004.....18	730144-003.....19	730150-002.....19
730058-002.....17	730066-008.....18	730527-005.....18	730144-004.....19	730150-003.....19
730058-003.....17	730066-012.....18	730527-008.....18	730144-005.....19	730150-004.....19
730058-004.....17	730066-052.....18	730527-012.....18	730144-008.....19	730150-005.....19
730058-005.....17	730066-053.....18	730527-052.....18	730144-012.....19	730150-008.....19
730058-008.....17	730066-099.....18	730527-053.....18	730144-052.....19	730150-012.....19

TMS Part #.....Page #	TMS Part #.....Page #	TMS Part #.....Page #	TMS Part #.....Page #	TMS Part #.....Page #
730150-052.....19	730156-012.....19	730147-008.....20	750053-225.....21	750133-400.....23
730150-053.....19	730156-052.....19	730147-012.....20	750054-225.....21	750133-500.....23
730150-099.....19	730156-053.....19	730147-052.....20	750055-125.....21	750133-600.....23
730151-002.....19	730156-099.....19	730147-053.....20	750056-125.....21	750133-700.....23
730151-003.....19	730350-002.....19	730147-099.....20	750057-225.....21	750133-900.....23
730151-004.....19	730350-003.....19	730148-002.....20	750058-225.....21	750134-200.....23
730151-005.....19	730350-004.....19	730148-003.....20	750059-110.....21	750134-300.....23
730151-008.....19	730350-005.....19	730148-004.....20	750060-110.....21	750134-400.....23
730151-012.....19	730350-008.....19	730148-005.....20	750061-210.....21	750134-500.....23
730151-052.....19	730350-012.....19	730148-008.....20	750062-210.....21	750134-600.....23
730151-053.....19	730350-052.....19	730148-012.....20	750063-110.....21	750134-900.....23
730151-099.....19	730350-053.....19	730148-052.....20	750064-110.....21	750136-200.....23
730445-002.....19	730350-099.....19	730148-053.....20	750065-210.....21	750136-300.....23
730445-003.....19	730351-002.....19	730148-099.....20	750066-210.....21	750136-400.....23
730445-004.....19	730351-003.....19	730149-002.....20	750067-125.....21	750136-500.....23
730445-005.....19	730351-004.....19	730149-003.....20	750068-125.....21	750136-600.....23
730445-008.....19	730351-005.....19	730149-004.....20	750069-225.....21	750136-900.....23
730445-012.....19	730351-008.....19	730149-005.....20	750070-225.....21	750137-100.....23
730445-052.....19	730351-012.....19	730149-008.....20	750071-125.....21	750137-500.....23
730445-053.....19	730351-052.....19	730149-012.....20	750072-125.....21	750137-900.....23
730445-099.....19	730351-053.....19	730149-052.....20	750073-225.....21	750138-100.....23
730152-002.....19	730351-099.....19	730149-053.....20	750074-225.....21	750138-500.....23
730152-003.....19	730446-002.....19	730149-099.....20	750512-110.....21	750138-900.....23
730152-004.....19	730446-003.....19	730412-002.....20	750513-210.....21	750151-250.....23
730152-005.....19	730446-004.....19	730412-003.....20	750514-110.....21	750151-300.....23
730152-008.....19	730446-005.....19	730412-004.....20	750515-210.....21	750151-350.....23
730152-012.....19	730446-008.....19	730412-005.....20	750516-125.....21	750289-200.....24
730152-052.....19	730446-012.....19	730412-008.....20	750517-225.....21	750289-500.....24
730152-053.....19	730446-052.....19	730412-012.....20	750518-125.....21	750289-600.....24
730152-099.....19	730446-053.....19	730412-052.....20	750519-225.....21	750289-900.....24
730153-002.....19	730446-099.....19	730412-053.....20	750520-110.....21	750290-200.....24
730153-003.....19	730447-002.....19	730412-099.....20	750521-210.....21	750290-500.....24
730153-004.....19	730447-003.....19	730413-002.....20	750522-110.....21	750290-600.....24
730153-005.....19	730447-004.....19	730413-003.....20	750523-210.....21	750290-700.....24
730153-008.....19	730447-005.....19	730413-004.....20	750524-125.....21	750290-900.....24
730153-012.....19	730447-008.....19	730413-005.....20	750525-225.....21	750399-200.....24
730153-052.....19	730447-012.....19	730413-008.....20	750526-125.....21	750399-500.....24
730153-053.....19	730447-052.....19	730413-012.....20	750527-225.....21	750399-600.....24
730153-099.....19	730447-053.....19	730413-052.....20	750285-001.....22	750399-900.....24
730154-002.....19	730447-099.....19	730413-053.....20	750285-002.....22	750400-200.....24
70154-003.....19	730448-002.....19	730413-099.....20	750286-001.....22	750400-500.....24
730154-004.....19	730448-003.....19	730483-001.....20	750286-002.....22	750400-600.....24
730154-005.....19	730448-004.....19	730484-001.....20	790078-001.....22	750400-900.....24
730154-008.....19	730448-005.....19	730489-002.....20	790078-002.....22	750401-200.....24
730154-012.....19	730448-008.....19	730489-003.....20	790078-003.....22	750401-500.....24
730154-052.....19	730448-012.....19	730489-004.....20	790078-004.....22	750401-600.....24
730154-053.....19	730448-052.....19	730489-005.....20	790078-006.....22	750401-900.....24
730154-099.....19	730448-053.....19	730489-008.....20	790079-001.....22	750402-200.....24
730155-002.....19	730448-099.....19	730489-012.....20	790079-002.....22	750402-500.....24
730155-003.....19	730449-002.....19	730489-052.....20	790079-003.....22	750402-600.....24
730155-004.....19	730449-003.....19	730489-053.....20	790079-004.....22	750402-900.....24
730155-005.....19	730449-004.....19	730489-099.....20	790079-005.....22	770001-001.....24
730155-008.....19	730449-005.....19	750043-110.....21	790079-006.....22	770003-003.....24
730155-012.....19	730449-008.....19	750044-110.....21	790079-007.....22	770003-007.....24
730155-052.....19	730449-012.....19	750045-210.....21	790079-008.....22	770002-001.....25
730155-053.....19	730449-052.....19	750046-210.....21	790079-010.....22	770080-001.....25
730155-099.....19	730449-053.....19	750047-110.....21	790122-002.....22	780057-002.....25
730156-002.....19	730449-099.....19	750048-110.....21	790122-003.....22	780060-001.....25
730156-003.....19	730147-002.....20	750049-210.....21	790122-004.....22	730286-001.....18/19
730156-004.....19	730147-003.....20	750050-210.....21	790122-006.....22	730286-002.....18/19
730156-005.....19	730147-004.....20	750051-125.....21	750133-200.....23	720750-001.....7/8
730156-008.....19	730147-005.....20	750052-125.....21	750133-300.....23	



2020 FAX ORDER FORM

FOR YOUR CONVENIENCE, THIS FORM MAY BE COPIED & RE-USED

BILL TO:

Customer Account Number: _____

Customer Account Name: _____

Address: _____

City: _____ State: _____ Zip: _____

Phone: (____) _____

Fax: (____) _____

PO #: _____ Your Name: _____

SHIP TO: (IF DIFFERENT THAN BILL TO)

Customer Name: _____

Attention: _____

Department: _____

Street Address: _____

City: _____ State: _____ Zip: _____

PRODUCT SELECTION

QUANTITY	PRODUCT NUMBER	PRODUCT DESCRIPTION	PRICE EACH	TOTAL

METHOD OF PAYMENT

☐ Check (Payable To TruMedical Solutions, LLC)

☐ Credit Card

☐ Visa

☐ Mastercard

☐ Discover

☐ American Express

Account number

--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--

Expiration Date: _____

Name (Please Print): _____

Signature: _____

Date: _____

Credit Card Billing Address: _____

City: _____ State: _____ Zip: _____

Sales tax applies where required by law. If you are tax exempt, please fax your tax exemption certificate along with your order.

Thank You For Choosing TruMedical Solutions, LLC

MERCHANDISE TOTAL _____

SHIPPING & HANDLING _____

SUBTOTAL _____

SALES TAX (On Subtotal) _____

TOTAL _____



Terms & Conditions

- For your convenience, we accept Visa, MasterCard, Discover and American Express.
- Terms of sale are Net 30 days upon credit approval.
- Please have correct part numbers and customer account number ready when ordering.
- A 1.5% finance charge per month may be added to all past due invoices.
- All products, terms and prices are subject to change without prior notice. Prices in U.S. Dollars.
- While every effort has been made to ensure the accuracy of the information provided within this publication, TruMedical Solutions, LLC. reserves the right to correct any errors and/or omissions found herein.
- All prices are FOB shipping point and all shipping charges are the responsibility of the purchaser and will be added to the invoice.

Return Goods Policy

- All returns must have a RGA number (Return Goods Authorization). A return authorization number may be obtained through our Customer Service Department.
- Products returned to TruMedical Solutions, LLC. 30 days after shipment shall be subject to a 30% restocking charge.
- No returns, exchanges or refunds will be accepted on custom or special ordered items.
- All returned items must be in saleable condition, accompanied by the original packaging material.

Warranty

- TruMedical Solutions, LLC. honors the warranty offered by the original manufacturer.

Motor Freight Policy

- Merchandise is carefully packed and inspected before leaving our facility. Responsibility for its safe delivery is assumed by the carrier upon acceptance of the shipment. Claims for loss or damage sustained in transit must be made to the Motor Freight Carrier.
- Visible loss or damage must be noted on the freight bill and signed by the carrier's agent. Failure to adequately describe external evidence of damage may result in the carrier refusing to honor a damage claim.
- Content of package(s) must be inspected immediately
- Concealed loss or damage means loss or damage which does not become apparent until the merchandise has been unpacked. When damage is discovered, submit a written request for inspection by the carrier's agent within fifteen days of the delivery date. All packing materials must be saved for inspection so that a proper claim may be filed. Damage is the carrier's responsibility. By following the instructions, we guarantee our full support of your claims to protect you against loss or concealed damage.

ORDERING INFORMATION

Order Options:

Order by phone, fax, mail or web.

Phone Toll Free:

1-877-882-7844

Hour of Operation:

Monday thru Friday

8:30 am to 6:00 pm EST.

Fax Toll Free:

1-877-882-7845

24 hours a day, 7 days a week.

Mail To:


TruMedical Solutions, LLC.

PO Box 1869

Collegedale, TN 37315

Order By Web:

www.tru-medical.com



We Accept All Major Credit Cards:





TRUMEDICAL SOLUTIONS IS A FULL LINE DISTRIBUTOR FOR:

3-Point Products

Altimate Medical, Inc.
Andover Healthcare
Armedica Manufacturing Corporation
Axelgaard Manufacturing Co., Ltd.
Bailey Manufacturing
BioMedical Life Systems
Biotone
BSN Medical Inc.
BV Medical
Chesapeake Medical
Compass Health
Contemporary Design Company
Core Products International, LLC
Cramer Products Inc.
Defibtech, LLC
DeRoyal Industries
DJO, LLC.
Dycem USA
Dynatronix Corporation
Earthlite Massage Tables

Fabrication Enterprises

FabriFoam Products
Fitness Systems, LLC.
Fitter International
GPI Anatomicals
Graham Professional Medical
Harbor Linen
Hausmann Industries, Inc.
Ideal Medical Products, Inc.
Lafayette Instrument Company, USA
Lohmann & Rauscher, Inc.
Magister Corporation
Matrix Fitness
Mediflow, LLC.
Mettler Electronics Corp.
N Z Manufacturing, Inc.
NK Products
Nonin Medical, Inc.
North Coast Medical, Inc.
OPTP
Ossur Americas

Pivotal Health Solutions

PrePak Products, Inc.
ProBand Sports Industries, Inc.
RehabPro, Inc.
Rich-Mar, Inc.
Safetec of America, Inc.
SCIFIT
Silipos, Inc.
Spirit Fitness
SportsArt America, Inc.
Squid Compression
Swede-O Inc.
The Hygenic Corporation
The Joint-Jack Company
The Prometheus Group
Total Gym Commercial, LLC.
Travanti Medical
USA Sports, Inc.
Whitehall Manufacturing
Wolters Kluwer Health

TruMedical offers a full line of capital equipment & supplies to meet the need of every clinic & clinician.

To view a complete product listing, go to our website at www.tru-medical.com.

If you have any questions call our customer service at 877-882-7844.



PO Box 1869 • Collegedale, TN 37315 • P: 877-882-7844 • F: 877-882-7845 • www.tru-medical.com