



## 3-Section Hi-Low Treatment Table



The Tru™ 3-Section Hi-Low Treatment Table is built to last. The frame is built of high quality square section steel tubing for maximum rigidity and stability. The table is also equipped with standard nose hole (and nose hole plug) for patient comfort while in the prone position. The Tru table offers an attractive closed base system that hides all electrical components making the table look spotless. This also makes your Tru table an absolute breeze to clean. The Tru 3-Section Hi-Low Treatment Table exhibits functionality & refinement.

### Item numbers:

**730459-xxx** Tru™ 3-Section Hi-Low Treatment Table w/ Footswitch control

**730460-xxx** Tru™ 3-Section Hi-Low Treatment Table w/ All-Round Bar control (ARB)

### Specs:

Overall Length: 77.75"L

Overall Width: 26"W

Head Section: 16.25"L x 26"W

Middle Section: 26"L x 26"W

Foot Section: 34.75"L x 26"W

Max Height: 37"

Min Height: 17.75"

Hi-Low Controls: Foot Switch or All-Round Bar (ARB)

2" Firm Density Foam Tops

Seamless Corners on Upholstered Components

Offered in (8) Upholstery Colors

Integrated, Retractable 3" Caters

110 Volts, 50-60Hz Motor

500 lb. Lifting Weight Capacity

Table Shipping Weight: 232 lbs.

Table Weight: 210 lbs.

Two Year Warranty